

Variant Analysis with UMI for Long-read Technology (VAULT)

[Edit](#)[Manage topics](#)

2 commits

1 branch


0 packages

0 releases

GPL-3.0

Branch: master ▾

[New pull request](#)[Create new file](#)[Upload files](#)[Find file](#)[Clone or download ▾](#)

 Chongwei add pic of 5' 3' umi		Latest commit 598bd45 7 minutes ago
example	add pic of 5' 3' umi	7 minutes ago
tools	Initial commit	13 minutes ago
.DS_Store	Initial commit	13 minutes ago
.gitattributes	Initial commit	13 minutes ago
LICENSE	Initial commit	13 minutes ago
README.md	add pic of 5' 3' umi	7 minutes ago
__init__.py	Initial commit	13 minutes ago
_version.py	Initial commit	13 minutes ago
check_umi.sh	Initial commit	13 minutes ago
variants_calling.py	Initial commit	13 minutes ago
vault.py	Initial commit	13 minutes ago

 README.md

VAULT

Variant Analysis with UMI for Long-read Technology (VAULT)

VAULT is a tool for analyzing UMI-labeled reads, works for both error-prone long reads and accurate single-end/paired-end short reads.

More detail: [Long-read Individual-molecule Sequencing Reveals CRISPR-induced Genetic Heterogeneity in Human ESCs](#)

INSTALL

Prerequisites

Anaconda (<https://www.anaconda.com/distribution/>) or Miniconda (<https://conda.io/miniconda.html>). For example, users may download and install the following Anaconda3 package:

```
wget https://repo.anaconda.com/archive/Anaconda3-2019.10-Linux-x86_64.sh
bash Anaconda3-2019.10-Linux-x86_64.sh
```

Download the VAULT package

```
git clone https://github.com/milesjor/VAULT.git
cd ./VAULT/
```

Install all required modules

```
conda env create --name vault --file ./tools/vault_env.yaml
conda activate vault
python -m pip install numpy
python -m pip install pandas
```

Parameters

The parameters of VAULT

Users can use the following command to print out the parameters:

```
python ./vault.py -h
```

A list of available parameters:

```
optional arguments:
  -h, --help            show this help message and exit

Required options:
  -u UMI_ADAPTER, --umi_adapter UMI_ADAPTER          # UMI primer sequence
                                                         UMI sequence, automatically detect \
                                                         NNNATGCNNN as UMI
  -s SAVE_PATH, --save_path SAVE_PATH                path/to/save/
  -r REFER, --refer REFER                            # the reference for alignment, \
                                                         should be the amplicon sequence
                                                         path/to/ref.fa
  -q FASTQ, --fastq FASTQ                            path/to/reads.fastq, or fastq.gz

Optional options:
  -e ERROR, --error ERROR                            # related to raw read error rate
                                                         error tolerance rate for umi analysis [0.11]
  -t THREAD, --thread THREAD                          # thread for processing UMI
groups                                                thread/process number [5]
  -T THRESHOLD, --threshold THRESHOLD                 # ignore UMI group with read
number <= [5]                                       Threshold of read number for snp analysis [5]
  -b BASH_THREAD, --bash_thread BASH_THREAD           # thread in variant calling
                                                         Thread for running bash cmd [1]
  -p PE_FASTQ, --pe_fastq PE_FASTQ                   read2.fastq for illumina pair-end sequencing
  -a {sr,map-ont,map-pb}, --align_mode {sr,map-ont,map-pb} # [sr] for Illumina, [map-ont]
for Nanopore,\                                     [map-pb] for PacBio
                                                         parameter in alignment, minimap2 -ax \
                                                         [sr|map-ont|map-pb]
  --unmapped_reads extract mapped reads before UMI analysis # will filter out reads that
cannot align to reference
  -v, --version      show the current version
```

Example

Command

```
python ./vault.py -u CATCTTACGATTACGCCAACCCTGCGGNNNNNTGNNNNNGACACATTCTCCCAGGCCCTACTT \
-q ./example/nanopore_reads.fastq.gz \
-s ./example/result \
-r ./example/reference.fa \
-e 0.11 \
-a map-ont \
--unmapped_reads \
-t 4 \
-b 1
```

Results

```
./result/
├─ nanopore_reads_alignment_summary.log # raw reads alignment summary
├─ nanopore_reads.bam
├─ nanopore_reads.bam.bai
├─ nanopore_reads_mapped.fastq # reads that can map to reference
├─ nanopore_reads_mapped.lst
├─ nanopore_reads.sam
├─ grouped_reads # fastq reads for every UMI groups
│   └─ perfect_umi
├─ snp
│   ├── all_snp_from_perfect_umi.vcf # raw variant calling result for small variants
│   ├── all_sv_from_perfect_umi.vcf # raw variant calling result for large variants (>=30bp)
│   └─ flt_sv_from_perfect_umi.vcf # filtered variant calling result for large variants
(>=30bp)
│   ├── pass_snp_from_perfect_umi.vcf # filtered variant calling result for small variants
│   ├── all_sv_from_perfect_umi.filtered.0.5.vcf # customized SV filter result
│   └─ all_sv_from_perfect_umi.filtered.0.5.sorted.vcf # customized SV filter result and sort by
position
│   └─ perfect_umi
├─ umi_analysis # UMI analysis result for 5' and 3' end of reads
│   ├── 3end_UMIs
│   └─ 5end_UMIs
```

Individual UMI group folder

```
./result/snp/perfect_umi/
├─ 14_ATCGATGATTTT_AAAATCATCGAT # 14 reads in this group, 5' UMI is ATCGATGATTTT, 3' is
AAAATCATCGAT
├─ 33_GACATTGTCTGG_CCAGACAATGTC
├─ 35_5end_AACAGTGCTGCT # 35 reads in this group, all reads from 5' UMI
├─ 5_3end_AAAAACATGGCA # 5 reads in this group, all reads from 3' UMI
├─ 7_5end_ATTCTTGGTGTC
├─ 7_CTATGTGAAGAA_TTCTTCACATAG
├─ 8_3end_ACAAGCAAAAAA
├─ 8_AGTTGTGCCATA_TATGGCACAAC
├─ 8_CCGCGTGAGATG_CATCTACGCGG
├─ 8_CGTTGTGTTACT_AGTACACAACG
├─ 8_CTATTTGTCACT_AGTGACAAATAG
├─ 8_GGGTTTGGTTTG_CAAACCAAACCC
├─ 8_GTGGGTGACGGG_CCCGTACCCAC
├─ 8_GTGTTTGTAGT_TCTAACAACAC
├─ 8_TCAATTGCAGAA_TTCTGCAATTGA
├─ 8_TTACTTGATTTT_AAAATCAAGTAA
├─ 8_TTGGATGGAAGT_ACTTCCATCCAA
├─ 9_AAAGATGCGCGT_ACGCGCATCTTT
├─ 9_AGAAATGATAGC_GCTATCATTCT
├─ 9_ATCGATGGTGCG_CGCACCATCGAT
├─ 9_ATGTTTGCCAAT_ATTGGCAAACAT
├─ 9_CTAAGTGCTTAT_ATAAGCAGTTAG
├─ 9_GAGAATGAGTAC_GTACTCATTCTC
├─ 9_GTTTATGTACAT_ATGTACATAAAC
```

5' end UMI and 3' end UMI

